

Standard Operating Procedure (COVID-19)

AREA	DEVICE	LEVEL OF DISINFECTION	FREQUENCY	WHO	DISINFECTING AGENT
Exam Rooms	Spuds, Alger Brush, Lacrimal Dilator, Cannulas	HIGH	1	a	B
	Tonometer/Pachimeter probes	HIGH	1	a	B
	Gonioscopy Lenses	HIGH	1	a	B
	Contact Lenses	HIGH	1	a	B
	Diagnostic Equipment (phoropter, perimeter, OCT, preliminary testing)	LOW	1	b	A
	Occluders / Eye patches	LOW	1	a	A
	Sinks	LOW	3	b	A
	Exam Chair and Unit	LOW	1	a	A
	R/G Glasses	LOW	1	a	A
	Trial Frame	LOW	1	a	A
	Hand Held Instruments	LOW	1	a	A
Laboratory	Frame Warmer	LOW	2	b	A
	Frames on Display	LOW	1	b	C
	Lab Hand Tools	LOW	2	b	A
Administrative	All countertops	LOW	2	b	A
	Computer keyboards, mouse and telephones	LOW	3	b	A
	Credit Card Terminal	LOW	1	b	A
	Staplers	LOW	3	b	A
	Pens	LOW	3	b	A
	Printer	LOW	3	b	A
	Waiting Area	LOW	1	b	A
	Door Handles	LOW	3	b	A
	Washrooms	LOW	3	b	A
Light Switches	LOW	3	b	A	

Frequency Legend	Disinfecting Agent	Who Legend
1. Before/After Direct Patient Contact	A. CaviWipes / 70% Alcohol	a. Optometrist
2. Numerous Times a Day	B. 6% Peroxide	b. Staff
3. End of Every Day	C. Soap/Water	
4. Weekly		
5. Monthly		